English

**DEX**cam™ 4HD

# Intra-oral Camera

User Manual

# DEXcam4 HD Intra-oral Camera Features

The DEXIS DEXcam4 HD intra-oral camera provides live high definition video display, as well as freeze and capture capabilities, for dental imaging. The camera is not meant to be used for diagnostic purposes. There are no contraindications for use of the camera.

- Proprietary, seven-element glass optical lens produce exceptional image quality
- Uses an HD Sony CMOS sensor with resolution of 1.3 Mega-pixels
- Two button capture control capability directly on the handpiece
- One-touch focus adjustment
- Detachable cable for easy sharing
- Robust aluminum construction.

## Intended Use

This intra-oral camera system and accessories are indicated for use to provide the dentist and the patient with a view of the mouth before and after the dental procedure, which assists the dentist in describing the dental procedure being performed as well as showing the results.

- This camera is designed to be operated by healthcare professionals who are educated and competent in the techniques described in this manual.
- Visual impairment is permissible so long as the user can resolve the necessary details of the accompanying documentation, user interface and image data. Use of corrective measures such as glasses or contact lenses is permissible.
- Audio impairment is permissible as there are no sounds emitted from the device.
- General physical impairments involving the arms, legs and/or motor skills are permissible so long as the user can perform all of the tasks required for proper operation of the device as described in this Manual.

## CONVENTIONS USED IN THE MANUAL

The following conventions are used to bring the operator's attention to important information:

 <b>WARNING</b>	Alerts the operator that failure to follow the procedure could cause serious bodily injury or death.
 <b>CAUTION</b>	Alerts the operator that failure to follow instructions may result in minor or moderate bodily injury or damage to the device.
<b>Please Note:</b>	Highlights important or unusual points.

## Proficiency with your New Device

Please become acquainted with your new DEXcam4 HD intra-oral camera, it is your tool. As with all new clinical tools, it is necessary to invest a certain amount of time for practice in order to become proficient. Learn to use it well and it will become an effective aid.

Please note that this manual assumes that new users possess basic computer skills and an understanding of the Windows® operating system. Absent this experience, we highly recommend that you obtain these skills through a computer course, video, or textbook. Your DEXIS representative may be able to suggest (although not endorse) one or more of these computer learning resources.

## DEXcam4 HD Intra-oral Camera Specifications


Name	DEXcam4 HD
Power	DC 5V, 500mA
Image Sensor	1/4" SXGA EDR 1.3 Mega-pixels
Resolution	1280 (H) x 720 (V)
AGC/ White balance	Auto
Lamp	High luminance white chip LED: 4EA
Cable	Length: 8.9 ft. (2.7 m) Connector: USB

## Package Contents

## Names and Functions of the Camera Parts

This device is an intra-oral camera system that provides the dentist and the patient with intra-oral and extra-oral images. This assists the dentist in describing the dental procedure being performed, as well as allowing the dentist to show the results to the patient.

	<p>1</p> <p>Viewing Head</p>	<p>Seven-element glass optical lenses</p>
	<p>2</p> <p>One-touch Focus Adjustment</p>	<p>5 position zoom stepper switch offers a focus mechanism to view intra-oral and extra-oral images without a lens change. The face icon captures smile, arch, and portrait images. The tooth icon captures up to three teeth. The magnifying glass captures a single tooth, detailed image,</p>
	<p>3</p> <p>Dual Capture Buttons</p>	<p>Freeze live images with a click of the two capture buttons.</p> <p>- Single-touch: Freeze image</p>
	<p>4</p> <p>Light Source</p>	<p>Four bright white LEDs</p>
	<p>5</p> <p>Quick-disconnect Cable</p>	<p>Easy sharing using multiple operatories.</p>

# DEXcam4 HD Installation

## Minimum system requirements:

- Microsoft Windows 7 or higher
- DirectX 9.0b or higher
- Pentium 4 at 2 GHz or higher
- 1 GB RAM
- Monitor with minimum 1024 x 768 resolution capability
- Graphics / display card capable of minimum 1024 x 768 resolution at 24-bit color and 32 MB video RAM
- CD-ROM drive
- USB 2.0 port

The installation procedure has two parts:

First - install DEXcapture software (refer to the DEXcapture Software manual for the software installation).

Second - install the intra-oral camera.

**Important!** Be sure to follow the instructions to install the software before connecting the camera. Otherwise the camera will not function properly.

## Install the Intra-oral Camera

### **WARNING**

No modification of this device is allowed.

### **WARNING**

Loading incompatible software or viruses may result in loss of patient data.


### **WARNING**

Do not connect multiple cameras to the same acquisition computer as this may lead to unexpected behavior by the software.

1. After installing DEXcapture, connect the USB cable of the camera to an available USB port on your computer.

**NOTE:** Remove the protective film covering the glass of the viewing head before using the camera.

2. Launch DEXcapture module from DEXIS software by clicking on the DEXcam icon.
3. Click the Display Settings slider button to display Application and Device settings.
4. Click the DEXcam4 HD device icon to connect the camera.
5. Refer to the DEXcapture Software Manual for all other software related settings.

6. Install the camera holster in a place that is convenient and does not exceed the maximum length of the USB cable, which is 8.9 ft. (2.7 m). The camera holster should also be positioned so that the cable from the camera does not create a tripping hazard. To install the holster, clean the mounting surface, remove the backing from the double sided-tape on the holster, then firmly press the holster against the mounting surface. After approximately two hours, the holster can be used to store the camera.

## Intra-oral Camera Sheath

Refer to the intra-oral camera sheath instructions on the sheath box for installation and replacement procedures.

To order more sheaths, please contact your local sales representative.

### **WARNING**

To help prevent cross-contamination between patients, place a new hygienic sheath on the camera for each new patient.

### **CAUTION**

The sheath must be removed after each patient. The intra-oral camera must be cleaned and disinfected after each patient. After cleaning the camera and handpiece, install a new camera sheath.

Intra-oral camera sheaths are applied parts. After each replacement, check the camera sheath for any tears or rips. Replace any damaged sheaths before use.

Avoid contact between the camera sheath and any sharp objects, such as the patient's teeth, orthodontia, etc. If there is any form of contact made, inspect the sheath for any tears or rips and replace it if necessary.

The camera sheath is disposable and should be replaced before each use of the camera. Dispose of sheaths following the normal dental office procedure for biomedical waste.

The camera sheath is designed only for intra-oral use. If you use it for taking extra-oral images, the images will not be clear as compared to intra-oral use because the transparent film reflects light and changes the focus.

# DEXcam4 HD Operation

To use the camera, remove it from the holster. It turns on and begins to send video images.

1. Preset the camera to view either an intra-oral or extra-oral image.

## CAUTION

Be sure to use a new sheath when capturing intra-oral images.

2. Position the camera.
3. Adjust the camera focus.
  - Press + to zoom in, which moves the LED indicator to the right.
  - Press - to zoom out, which moves the LED indicator to the left.
4. Press and release the two Capture buttons to freeze the image on your viewing screen.

## CAUTION

Ensure that the camera is properly placed into the holster with the cord end facing down. The camera will fail to turn OFF if placed in the holster with the cord end facing up.

5. Place the camera in the camera holster to turn the camera off.

## Handling

After extensive use, the camera focusing motor may need to be reset (only as needed).

1. Press the - button to move the LED indicator to the leftmost position.
2. Then press the + and - buttons simultaneously until five LED indicator lamps are lit then turned off.
3. Using the + and - buttons, focus on an object more than 10 feet away.
  - Press + to zoom in.
  - Press - to zoom out.
4. Once the image is focused and clear, press the + and - buttons simultaneously until the five LED lamps blink again.
5. The Camera focus is now reset.

## CAUTION

Failure to follow instructions may result in damage to the product.

- Do not drop or strike the camera since the internal precision optics could be damaged.
- When the camera is not in use, place the camera in the camera holster, which will turn it off.
- The intra-oral camera system should only be used for dental purposes.

- The camera is not waterproof. Do not immerse in water and avoid situations where water or other foreign matter could get on or inside the product. Failure to do so may damage the product.
- If the camera is wet, or if your hands are wet, do not use the product.
- Avoid subjecting the camera to vibrations or shocks. When the product is not in use, make sure it is placed in the camera holster.
- Examine the camera housing for visible damage. Be careful not to use the camera if it is damaged. Doing so may cause injury.
- Avoid using the camera in a dusty environment.
- Remove the protective film covering the glass of viewing head.

## Important Safety Precautions

### CAUTION

Failure to follow instructions may result in damage to the product.

- If the camera gets wet, disconnect it, wipe it down, then allow it to air dry. Do not attempt to dry with a heater or microwave.
- During severe thunderstorms, stop using the camera. Failure to do so may cause damage to the product.
- Do not disassemble, drop, or otherwise subject the camera to hard shocks. It may cause electric shock, short circuit, or fire. The warranty is null and void in such cases.

Attention: US federal law restricts the sale of this device to licensed dental practitioners only.

## Cleaning and Disinfecting

### Important:

Before applying cleaning or disinfecting solution:

- disconnect the camera from power, and
- disconnect the USB cable from the USB port.

Do not soak the camera.

Always dry completely after cleaning or disinfecting.

To clean: Use a lint-free cloth to apply a solution of mild soap and water.

To disinfect: Use a lint-free cloth to apply disinfectant and rub for at least on minute, using one of the following disinfectants:

- Isopropyl alcohol (70%)
- CaviCide
- CaviCide1
- CaviWipes
- CaviWipes1

**NOTE:** After applying the disinfectant to the cotton cloth, excess solution must be squeezed out before wiping the camera. When using the wipes, excess solution must also be squeezed out before wiping the camera.

## Label Location

The following label is located on the USB cable of the camera.

## Product Symbols

Indicates the product should be used only once. The symbol is found on the packaging for the sheaths.

This symbol on the product and/or accompanying documents means that used electrical and electronic products should not be mixed with general household waste.

Please refer to the written instructions of this manual.

Direct Current (DC)

The ETL listed mark indicates that Intertek has certified the product described herein under control number 4005847 to be in compliance with the applicable regulations. Intertek is:

- a Nationally Recognized Testing Laboratory by the Occupational Safety and Health Administration (OSHA) in the United States.
- a Certifying Body in Canada by the Standards Council of Canada.

The CE symbol represents that the product herein specified meets the provisions of European Council Directive 93/42 EEC concerning medical devices and Directive 2011/65/EU of the European Parliament and of the Council on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS).

Date of Manufacture

Manufactured for

Authorized Representative in Europe

## Service and Proper Disposal

The service life of the camera is five years. There are no serviceable parts or maintenance on the device. Consult the device warranty for additional information.

Properly dispose of the camera at the end of its useful life. Do not mix with general household waste.

## Specifications

### **USB Compatibility**

USB 2.0 or later

### **Electrical Ratings**

5VDC, 500mA

### **Requirements for safety approvals**

This device must be connected to a personal computer that meets the standard specified in IEC 60950, revision 2003.

### **Ambient conditions:**

Transportation and storage temperature: -40°C to +70°C (-40°F to +158°F)

Operating temperature: +10°C to +40°C (+50°F to +104°F)

Relative humidity: 30% to 75%, non-condensing

Air Pressure: 700 to 1060 hPa

Maximum surface temperature of camera at maximum ambient: 59°C (138°F)

# Manufacturer’s Declaration

The intraoral camera emits electromagnetic energy at levels compliant to international standards; however, interference between the camera and other devices is possible if located in close proximity. The customer or the user of the camera can help prevent electromagnetic interference by locating sensitive devices away from the camera.

This camera, classified as Medical Electrical Equipment, requires special precautions regarding EMC and must be installed and put into service according to the EMC information provided in the accompanying product documentation. Portable and mobile RF communications equipment can effect Medical Electrical Equipment. This camera complies with EMC requirements when used with the cables and accessories supplied with the product. The use of accessories and cables other than those sold by DEXIS, LLC and specified as replacement parts for internal components, may result in increased emissions or decreased immunity of this camera. This camera should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, this camera should be observed to verify normal operation in the configuration in which it will be used.

Guidance and manufacturer’s declaration - electromagnetic emissions for all EQUIPMENT and SYSTEMS

Guidance and manufacturer’s declaration – electromagnetic emissions		
This camera is intended for use in the electromagnetic environment specified below. The customer or the user of this camera should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	This camera uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The camera is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded:
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	<b>Warning:</b> This camera is intended for use by healthcare professionals only. This camera may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the camera or shielding the location.

Guidance and manufacturer's declaration - electromagnetic immunity for all EQUIPMENT and SYSTEMS

<b>Guidance and manufacturer's declaration - electromagnetic immunity</b>			
This camera is intended for use in the electromagnetic environment specified below. The customer or the user of this camera should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD)	±6 kV contact	±2, 4 & 6 kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
IEC 61000-4-2	±8 kV air	±2, 4 & 8 kV air	
Electrical fast transient/burst	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
IEC 61000-4-4	±1 kV for input/output lines	±1 kV for input/output lines	
Surge	±1 kV line(s) to line(s)	±0.5 & 1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
IEC 61000-4-5	±2 kV line(s) to earth	±0.5, 1 & 2 kV common mode	
Voltage dips, short interruptions and voltage variations on power supply input lines	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0,5 cycle  40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles  70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles  <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 sec	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0,5 cycle  40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles  70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles  <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of this camera requires continued operation during power mains interruptions, it is recommended that this camera be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
IEC 61000-4-8			
NOTE $U_T$ is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - electromagnetic immunity for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration – electromagnetic immunity			
This camera is intended for use in the electromagnetic environment specified below. The customer or the user of this camera should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	<p>Portable and mobile RF communications equipment should be used no closer to any part of this camera, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P} \quad 80 \text{ MHz}-800 \text{ MHz}$ $d = 2,3\sqrt{P} \quad 800 \text{ MHz}-2.5 \text{ GHz}$ <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people</p>			
<p><sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this camera is used exceeds the applicable RF compliance level above, this camera should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this camera.</p> <p><sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and this camera			
This camera is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this camera can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this camera as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

## Accessories

Part	Description	Part Number
	Camera Holster	PLU7611H
	USB Quick-disconnect Cable	PLU7631C
	Intra-oral Camera Sheaths Case (5 boxes - 500 count)	PLU7611S

# Additional Help and Support

## DEXIS on the Internet

[www.dexis.com](http://www.dexis.com)

The DEXIS device described herein does not contain parts which can be serviced.



## DEXIS Customer Care Center (United States and Canada)

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Kaltenbach & Voigt GmbH