
Turbo

- ① EN Class B Sterilizer
- ② ES Esterilizador Clase B
- ③ FR Stérilisateur de classe B

INSTRUCTIONS FOR USE
INSTRUCCIONES DE USO
DIRECTIVES D'UTILISATION

INDICATIONS FOR USE

The Maxima Turbo is an air-removal (pre-vacuum) table-top steam sterilizer intended for use by a health care provider to sterilize medical products by means of pressurized steam. It is suitable for the sterilization of dental and medical instruments that are validated to be sterilized by steam. The Maxima Turbo has not been designed to sterilize liquid loads, bio-medical waste or materials not compatible with steam sterilization. The processing of such loads may result in incomplete sterilization and/or damage to the autoclave. Please refer to the table below for program name, load description, sterilization temperature, exposure time, drying time and maximum load.

MAXIMA TURBO

Program name	Load description	Sterilization temperature	Sterilization time	Drying time maximum load	Maximum load
134°C	solid objects, small porous objects, simple objects recessed, narrow-clearance items, dental handpieces, and textiles; wrapped and unwrapped	134°C (273°F)	4:00 minutes	3 minutes	0.5 Kg /1.1 lbs
121°C	solid objects, small porous objects, simple objects recessed, narrow-clearance items, dental handpieces, textiles, and plastics; wrapped and unwrapped	121°C (250°F)	30 minutes	5:30 minutes	0.5 Kg /1.1 lbs
134°C FAST*	solid objects, non-porous objects, simple instruments (such as scissors, handles, pliers, chisels, probes, etc.), and dental handpieces; unwrapped	134°C (273°F)	4:00 minutes	N/A	0.5 Kg /1.1 lbs

* - Immediate Use Steam Sterilization cycle

The sterilizer is suitable for use in the vicinity of other powered medical products.

The Maxima Turbo device may not be used to sterilize liquids, biomedical waste or pharmaceutical products.

The device is intended for professional use by properly trained staff only.

SYMBOLS USED ON THE DEVICE

This symbol is located on the front of the appliance, on top of the drawer front, and advises extreme caution due to the high temperature of the compartment and its immediate surroundings.

This symbol is located on the unit's rating plate and identifies its individual serial number.

This symbol is located on the unit's rating plate and identifies the year the unit was manufactured.

This symbol is located on the equipment rating plate and identifies the equipment manufacturer.

This symbol can be found on the unit's rating plate and requires the user to read and follow the information in this manual.

Used electrical and electronic equipment must not be placed, thrown away or stored with other waste. Used equipment should be delivered to a local collection point for used electrical equipment, which is registered with the respective environmental protection office and conducts selective waste collection.

INDICATIONS FOR USE

- Maxima Turbo complies with IEC 60601-1-2:2014 (Edition 4.0)
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally. Max power cord length is 160 cm.
- Use of accessories and cables other than those specified or provided by Maxima Turbo could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Maxima Turbo S, including cables specified by Maxima Turbo. Otherwise, degradation of the performance of this equipment could result.
- The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
- The user is responsible for the installation, correct operation and maintenance of the device in accordance with instructions provided in this user manual. If needed, contact the service or the supplier of the device.
- The sterilizer is not intended for sterilizing liquids, biomedical waste or pharmaceutical products.
- The sterilizer must not be used if explosive gases or vapours are present in the air.
- After the cycle is completed, the load is hot. Remove tools or packs from the chamber using appropriate thermal gloves or equipment that prevents burns.
- Do not remove the rating plate or any other elements of labeling from the device.
- Follow guidelines for preparing tools for sterilization.
- Pouring water or other liquids on the device may cause a short-circuit.
- Prior to inspection, maintenance or servicing, turn off the device and disconnect it from the power source.
- Servicing may only be performed by trained service personnel and using original replacement parts.
- The device is intended for indoor use only.
- Max working high elevation of the device is 2,000 m above sea level.
- Pollution Degree 2: Normally only nonconductive pollution occurs. Temporary conductivity caused by condensation is to be expected.
- Overvoltage Category II.
- If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.
- Use only detachable cord that rated equal or greater to the equipment electrical ratings.


Read this Operating Manual carefully before using this device. Install and operate the device strictly as specified herein. Comply with all safety requirements for the device. This will ensure proper and safe operation of this device. Any other application, inconsistent with this manual, may lead to dangerous accidents. Restrict unauthorised personnel access to the device and train the personnel handling the device. An operator of this device is any person who, by training, experience and knowledge of applicable reference standards, manuals and occupational health and safety regulations has been authorised for the essential operation with the device and who is capable of identifying and avoiding the hazards related to operation of this product.



Always append this Operating Manual with the device if transferred to a new owner. The Operating Manual contains detailed information about assembly, installation, initial start-up, use, repairs and maintenance of the device. If the device is used as intended, this Manual will provide sufficient guidance to qualified personnel. Keep this Operating Manual close to the device and easily accessible at all times. As required by continuous improvement of the product, the manufacturer has the right to amend this Manual or make changes to the device without prior notice. Henry Schein shall not be liable for damage incurred during the wait for warranty service, any damage to the Customer's property other than this device, or errors caused by improper installation and/or improper operation of the device.

Detailed recommendations, counterindications and warnings are described in the relevant sections.

PRODUCT STRUCTURE AND ACCESSORIES




DEVICE SETUP

Please read these operating instructions carefully before using the Maxima Turbo. During operation, observe the guidelines in these instructions as well as the applicable local health and safety rules. The Maxima Turbo unit is designed for self-assembly by the user and do not require specialized installation on site. The user is responsible for proper installation of the unit on site in accordance with these instructions.

HEPA filter installation.

To protect the component during transport, a HEPA filter is not installed in the device. It was placed inside the chamber, on a tray. Take it out of the chamber, take it out of the bag and by yourself, carefully screw it into the slot provided for this purpose (see fig.). The filter should be screwed in by hand (approx. 7 turns) until resistance is felt - further screwing in will damage the filter and immobilize the sterilizer.

Feed water quality

It is unacceptable to use any liquids, solutions, chemicals, additives to the feed water - the device may be powered only with demineralized or distilled water of appropriate quality.

No chemicals or additives must be added to the sterilization water, even if they are specifically intended for use in steam generators, steam production or for use as additives in sterilization, disinfection, cleaning or corrosion protection.

The use of incorrect quality water or solutions other than water reduces the effectiveness of sterilization and leads to damage to the device and loss of warranty.

The total mineral content of the sterilization water must be less than 10 ppm or, in the case of conductivity measurement, it must be less than 15 $\mu\text{S}/\text{cm}$ - therefore tap water cannot be used as feed water for the device.

The table below shows recommended water hardness and conductivity parameters. Permissible parameters of water used for sterilization

- Hardness $<0.02 \text{ mmol / l}$
- Conductivity (at 20°C) $<15 \mu\text{S / cm}$
- PH - from 5 to 7.5

- Chemical additives - none

Water Filter & Water Quality Indicator

To ensure the proper feed water quality a water filter and water quality indicator has been attached to the inlet line. Any usage of improper water will begin to discolor the indicator. Check the indicator frequently to ensure proper water is being used. The use of incorrect quality water or solutions other than water reduces the effectiveness of sterilization and leads to damage to the device and a loss of warranty.

If the indicator is completely discolored the unit will need to be serviced and the inlet line will need to be replaced.


Connecting the WATER IN feed water

- Connect the blue connection hose to the blue feedwater connection, which is marked on the rear panel of the WATER IN device. The tubing should be inserted into the connector, to the depth of the black line (placed on the tubing).
- The other end of the blue tube should be immersed in the tank with the feed water. The device is equipped with a water suction pump, so it is not necessary to place the water tank above or on the same level as the device.
- To secure and immobilize the water supply hose to the machine, use the plug supplied with the hose and insert the plug into the opening of the water supply tank. An unsecured tubing may jump out of the water and cause sterilization errors.
- To prevent the hose from coiling in the water tank, install the included weight with a rubber ring at 2 cm from the end of the hose (fig.)
- The minimum water load in the tank is 300 ml.
- Remember and check that the blue tube is always immersed in the water.
- Check the water level in the tank on a regular basis, depending on the frequency of the processes

WATER OUT Wastewater connection

- Connect the red wastewater hose to the orange wastewater connector that is marked on the rear panel of the WATER OUT device. The tubing should be inserted into the connector, to the depth of the black line (placed on the tubing).
- Wastewater should be discharged directly to the water and sewage system or to a special wastewater tank. If a tank is used, the other end of the red tube should be placed in the wastewater tank. To secure and immobilize the water drain hose from the machine, use the plug provided with the hose and insert the plug into the opening of the wastewater tank. The unsecured tubing may jump out of the tank and flood the room.
- Make sure that the red tube is never immersed in water, otherwise the water will not drain properly, causing sterilization errors.
- The wastewater tank or drainage outlet must always be located below the unit.
- When using wastewater tanks, we recommend that you use tanks with the same capacity as the feed water tank. Their simultaneous replacement will protect against the possibility of wastewater overflow.
- The water level in the tank should be checked on a regular basis, depending on the frequency of the processes.

Setting up the device


- The unit should be placed on a flat, level surface. Do not use the unit if it is at an angle.
- Only demineralized or distilled water should be used in the devices. The use of incorrectly specified water (see “Supply water quality”) reduces the sterilization efficiency and leads to loss of warranty due to damage to the device.
- Leave a space of 5 cm (2 in) at the rear of the appliance and at least 1 cm (0.4 in) at the sides of the appliance from walls or other elements for proper ventilation and fresh air access.
- The unit should be located to provide easy access to the main switch located on the rear panel of the unit.
- Do not place the unit near a sink with water or other places where the unit could be flooded - danger of short circuit.
- The device should be installed in a well-ventilated room, far from heat sources and rooms where mixtures of gases or liquids or other dangerous factors may arise.

The following conditions must be ensured for safe and efficient operation of the device:



- Ambient temperature from + 5 ° C (41°F) to + 25 ° C(77°F),
- Relative humidity from 0 to 90%,
- Storage temperature from -20 ° C (-4°F) to + 60 ° C (140°F),
- 0 to 90% relative humidity.

Correct connection of the hoses to the feed and wastewater connections

Correct installation of the weight with the securing ring

Correct positioning of the hoses in the water tank

PREPARATION AND LOADING OF INSTRUMENTS

The devices are suitable for sterilization of loads for which steam sterilization is stipulated. The instruments must be sterilized only if they are clean and dry. Therefore they must be washed and disinfected before being loaded onto the tray in accordance with applicable regulations. Residues of used chemicals or other solid particles may render the sterilization process impossible or even damage the device. In addition, sterilization of instruments that have not been previously cleaned and disinfected is a biological hazard and can lead to damage (both to the instruments and to the sterilizer). For instruments that must be lubricated, use steam sterilization lubricants. Excess lubricant must also always be removed.

- The 134 ° C (273°F) FAST program is intended only for sterilization of unwrapped instruments.
- After sterilization with 134°C (273°F) FAST the instruments are moist - it is recommended to leave the door open for a few minutes to evaporate excess moisture.
- After sterilization of unpackaged loads, they are intended for immediate use.

Preparation of packages

It is recommended to use sterilization packaging that meets the requirements of EN ISO 11607-1: 2020 and EN 868-2: 2017. Appropriate packaging should:

- Ensure good penetration of the sterilizing agent into the inside of the package.
- Provide resistance to damage during the sterilization process.
- Ensure tight and durable sealing of contents and safe removal for use.
- Provide a barrier to microorganisms and unwanted substances and contaminants.
- Use disposable sterilization sleeves (intended for steam sterilization) as recommended by their manufacturer.
- Sleeves should only be filled to about 3/4 full to allow for proper sealing and minimize the risk of damage to the package.
- A distance of at least 30 mm should be maintained between the seals and the sterilized load.
- The sharp edges of the load should be protected to avoid damaging the packaging.
- The packaging material must not be placed too loosely or be tightly stretched so as not to affect pressure variations during sterilization.
- a label with information about the contents of the package, the code of the packer, date of sterilization and use-by date as well as parameters of the sterilization process should be placed on the package.

Arranging packages on the sterilizer tray (sterilization of packed items)


- Packets on the tray should be placed so that the paper side touches the paper side, because penetration of the sterilizing agent and air exchange can only take place through the paper.
- The packages should be placed on the tray in such a way as to eliminate the contact of the package with the door seal and the phase of the sterilization chamber - failure to comply with this rule may cause the chamber to become unsealed and the sterilization cycle incorrectly carried out
- The edges of the packs must not protrude beyond the sterilizer tray, as this would cause the chamber to leak and the sterilization cycle to fail
- If the sterilizer chamber is heavily loaded, the first packs must face the foil side towards the bottom of the tray (this ensures faster and more efficient drying of the packs).

Arranging tools on the sterilizer tray (sterilization of loads without packaging)

- In the case of sterilization without packages - the instruments must be placed in such a way that they do not have direct contact with each other, no element of them falls into the openings of the tray, and does not rest against the edge of the tray or protrude beyond its outline.
- Failure to follow the above recommendations may cause permanent and irreversible damage to the sterilization chamber phase, which will result in the lack of tightness of the sterilizer and loss of warranty.

Do not exceed the maximum load weight: 500 g for the Maxima Turbo device. In order to verify the correctness of the performed sterilization, it is recommended to place a sterilization chemical indicator in the chamber for each process, which discoloration during properly conducted sterilization.

Sample placing of sterilization packages.


Example of solid load.

OPERATING THE DEVICE

First run

Turn on the device using the main switch located on the rear panel of the device. Make sure that the supply and wastewater hoses are connected correctly, and that there is water in the supply tank and the wastewater tank is empty. Place the load on the tray, slide the drawer in and lock it, i.e. turn the knob clockwise. The sterilizer provides audible information about changes and the completion of the process.

Program selection

Depending on the type of load to be sterilized, the user must select the appropriate program - in accordance with the instructions of the sterilizer manufacturer and the load manufacturer's recommendations.


Maxima Turbo program	134°C (273°F) FAST	134°C (273°F)	121°C (250°F)
Type of load	Unwrapped instruments only	Wrapped or unwrapped instruments	Wrapped or unwrapped instruments
Process temperature	134°C (273°F)	134°C (273°F)	121°C (250°F)
Pre-vacuum number	1	3	3
Sterilization duration	4:00 min	4:00 min	30 min
Drying duration	-	3 min	5:30 min
Total process duration*	100g: 7 min	100g: 13 min	100g: 45 min
Class	S	B	B
Load Description	Solid Objects, non-porous objects, simple instruments (such as scissors, handles pliers, chisels, probes, etc.) and dental handpieces; unwrapped	Solid objects, small porous objects, simple objects recessed, narrow-clearance items, dental handpieces, and textiles; wrapped and unwrapped	Solid objects, small porous objects, simple objects recessed, narrow-clearance items, dental handpieces, textiles, and plastics; wrapped and unwrapped


After switching on the device, the start screen appears on the display. To go to the next screen, press the screen once with your finger (anywhere).


After switching on the device, the start screen appears on the display. To go to the next screen, press the screen once with your finger (anywhere).

On the next screen, we can go to sterilization, test the device or open the information menu. The following sterilization programs can be selected from the main menu: 121 ° C (250°F) , 134 ° C (273°F) , and 134 ° C (273°F) FAST.



When the chamber is opened, the DOOR OPEN symbol flashes.

When the chamber is closed by turning the locking knob clockwise, the START symbol appears on the display, which informs that the chamber is properly closed.


If the USB memory is not inserted in the device, the USB disk symbol is not displayed in the lower right corner of the screen.

A message appears that USB memory is missing. The program data will not be saved. You can continue working without saving data to the USB memory stick by pressing the YES box, or abort the work by selecting the NO box to place the stick in the port and start the program steps from the beginning.

It is recommended to use a USB drive with each of the processes. The data stored on it will allow you to print reports from selected processes.

If work continues or the START field has been selected, a symbolic pressure diagram of the entire process is displayed on the screen with the current program progress highlighted against a background and information on the next process steps in the upper left corner of the screen. During the program run, the screen displays the temperature of the selected sterilization program or the current temperature in the process chamber in the lower left corner, the current pressure in the chamber in the lower right corner, and the time remaining until the end of the process in the upper right corner of the screen. This is an estimated time, which may be increased due to the weight and type of the load.

During the course of the program, in place of the START field, the STOP field is displayed, which allows the user to stop the process at any time. In the upper left corner of the screen, the names of the individual stages of the program are displayed consecutively, e.g. chamber lock, heating of the working chamber.

If the process was successful, the display alternately shows information screens about the end of the process and the sterility of the load and about the possibility of opening the chamber of the device. By pressing the FINISH field, you go to the welcome screen.


WARNING! After the process, the chamber, tray and load are hot. Be especially careful and use protective gloves to remove the load or wait until it cools down. In the 134C FAST program, instruments are hot and wet after sterilization.


Test programs

The devices are equipped with special test programs to verify the correct operation of the sterilizer.


Maxima Turbo	Bowie & Dick / Helix	Vacuum test
Process temperature	134°C (273°F)	-
Number of pre-vacuums	3	1
Sterilization time	3.5 min	-
Drying time	3 min	-
Total process time	15 min	16 min

By pressing the Test field you go to the test programs menu. From this level, we can choose between the Vacuum test and Helix / B&D test programs. You select the appropriate program by pressing the required field on the display.

When the device working chamber is closed, the word DOOR OPEN changes to START and by pressing this field the selected test program is started.

If a USB flash drive has not been inserted into the unit, the USB symbol is not displayed in the lower right corner of the screen and a message indicating no USB flash drive is displayed. The test program data will not be saved. You can continue without writing data to the USB flash drive by pressing the YES box or abort by selecting the NO box to place the USB flash drive in the port and start the program from the beginning.

If work continues or the START field has been selected, a symbolic pressure diagram of the entire process is displayed on the screen with the current program progress highlighted against a background and information on the next process steps in the upper left corner of the screen. During the program run, the screen displays the temperature of the selected sterilization program or the current temperature in the process chamber in the lower left corner, the current pressure in the chamber in the lower right corner, and the time remaining until the end of the process in the upper right corner of the screen. This is an estimated time, which may be increased due to the weight and type of the load.


Vacuum test

The vacuum test should only be performed on a cold device before starting work. During the test, the device verifies:

- Vacuum pump capacity.
- Tightness of the pneumatic system.


After selecting a vacuum test program and starting it with the button, the vacuum test program run screen is displayed, with Information about loss of pressure in the working chamber and the duration of the test.

[Vacuum test PASSED]

When the test program ran successfully.

[Vacuum test FAILED]

When the test program did not run successfully.

After pressing the CONTINUE field, the welcome screen is displayed.

During the vacuum test the sterilizer chamber must be completely dry and cold. If not, the vacuum test may not be reliable even if the sterilizer is fully operational. When the test is complete a message indicating the result will appear on the display. If the result is negative check, clean or replace the seal, clean the front edge of the chamber and repeat the test. If the Vacuum test fails again, contact the supplier or manufacturer.

Bowie & Dick test


The Bowie & Dick test, also known as the steam penetration test, simulates a small, highly porous load. To perform the test it is necessary to have a special test packet and placed inside the chamber. The package is not an accessory of the device, the user should purchase it on his own.

This test evaluates the device's performance in sterilizing loads of porous objects:

- Pre-vacuum performance and steam penetration.
- Temperature and pressure of saturated steam reached for a certain period of time.

How to run the test:

- The test must be carried out with an empty chamber in accordance with EN 13060.
- Place the Bowie-Dick test package in the chamber in the center of the tray.

After selecting the Helix / B & D test program and starting it with the START button, the program sequence screen is displayed.


Information about the process parameters is displayed.

The Helix / B & D test program can be stopped at any time by pressing the STOP field, which is associated with incorrect termination of the test.

When the test program is complete, alternating screens are displayed:

[READY] Test complete / [READY] Chamber is safe to open. The process chamber of the sterilizer can be safely opened.


The control test should be removed.

WARNING! The package will be hot.

To correctly interpret the test, refer to the instructions provided by the test manufacturer.

Open the package and remove the indicator chemical from the center of the package.

Positive result

The chemical indicator turned a uniform dark color over the entire surface.

Any color change, uneven coloring of the test, indicates the presence of air during the test cycle due to a malfunction of the sterilizer. If the test result is abnormal, check the expiration date of the test pack and repeat the test.

Negative result

A bright field remained in the middle of the test because of the remaining air inside the tested device.


Information menu

The info menu is accessible by pressing the Info. button.

This level displays information about the unit type, serial number, number of processes performed and available USB memory space for saving process data, as well as the COUNTERS service menu - process counters for filter changes. It is also possible to change the date and time. To set the date or time touch the digits on the display. When a certain field is selected, it starts flashing and the arrows for changing the value up or down are displayed. This is how to set the date and time correctly. Pressing the number again confirms it and you can move on to changing the settings of the next parameter. In the same way, you can select the language by clicking on its abbreviation.

The button labeled B turns off and on the blue backlight in the depth of the screen. Clicking the LED button launches the backlight control menu on the sides of the device.


Clicking a button turns off the button sounds.

LED lighting

LED lighting has two modes:

- Free mode, where the user (by moving the sliders) freely sets the colors, intensity and brightness of the light to their preference.
- Continuous mode that indicates the stages of the entire sterilization process with colors.

Counters

The Maxima Turbo sterilizer records the number of processes performed. This lets you know when you need to replace consumable parts.

Counter section with the counters reset to zero. The number of processes performed is on the left, while on the right is the number at which the component should be replaced or a service inspection should be performed. After replacing the filter, the user can reset the values by pressing the RESET button.

During the first process (Vacuum, Helix, FAST, 134/273, 121/250), the unit will record the current date as the date of the first process (bottom line on the counter screen). The device will count processes and days since the first startup.

Displaying the counter values in yellow or red does not block the operation of the unit. However, exceeding the required replacement time may have a significant impact on the operation of the unit and the sterility of the load. For replacement of individual components, please contact the manufacturer or supplier.

No.	Name	Replacement frequency (cycles)	Yellow Early replacement call (cycles)	Red Latest replacement after (cycles)
1	HEPA filter and Door with gasket	1000	980	1000



Restart

The restart of the process is forced when the user stops the process by pressing the STOP field, in the case of a power outage or an error during the process, for example lack of supply water.



If the STOP field is selected, the following messages are displayed alternately:

- The user has stopped the process.
- Equalizing the pressure in the working chamber.
- Process is incorrect, which means that the load is not sterile.

When the pressure in the working chamber is equalized, the following messages appear alternately on the display. Now you can open the device freely. After opening the chamber, the screen appears.

By selecting the field, we can return to the welcome screen. In the event of an error, we must additionally enter the 4-digit security code 0000. Entering this code is tantamount to the operator's declaration that he is aware that the sterilization process has not been carried out properly and that the batch is non-sterile.

DEVICE MAINTENANCE

In order to ensure the correct operation of the device, the user is obliged to perform the following maintenance actions.

Cleaning the tray

Keeping the tray clean helps to ensure proper operation of the appliance. The tray and its technical condition are a good indicator of using the correct water. A stony, brown tray indicates the use of poor quality water. It is recommended to clean the inside of the tray on a weekly basis with a mild detergent that does not add chlorine and does not react with aluminum. After cleaning, the tray requires thorough rinsing with water. In order to clean the tray properly, it should be removed from the front of the device.

Lift the tray up gently and pull away from the front face. The mounting pins have notches into which the drawer fits. Before reinstalling the tray into the unit, drain the tray and slide it over the front face pins and press down gently to lock.

Cleaning the process chamber

Keeping the chamber clean helps to ensure proper operation of the device. It is recommended to clean the inside of the process chamber once a week with a mild detergent without the addition of chlorine. After cleaning, the chamber should be wiped dry with a soft cloth. Cleaning should be performed on a cold chamber.

Cleaning of external surfaces

Cleaning the external parts of the device should be done with a soft cloth moistened with water and a mild detergent (no chlorine added and not reacting with plastics, varnish coatings, aluminum). Strong detergents should not be used. The use of mild detergents to maintain the device does not affect the risk of toxic components coming into contact with device components.

Cleaning the gasket

It is recommended to clean the gasket each time after 100 processes. Use warm, clean water and microfiber to clean the gasket (microfiber with silver particles is allowed). Blunt and sharp cleaning tools are not allowed. Cleaning with chemicals is not allowed. Carry out cleaning on a cooled down device, after opening the drawer. Be careful not to bend the drawer. After cleaning, leave the device open until the gasket is dry. During this time, protect the device against damage.

Replacement of consumables

The sterilizer is equipped with a high performance sterilization system. A message on the screen informs the user when each component should be replaced. If the sterilizer is in regular operation, alternating replacement screens will appear after the welcome screen has been pressed. The replace-

ment screens are described in detail in section “Warning Messages and Error Codes”.

Cleaning the water container

In order to ensure the correct parameters of the water supplying the device, it is recommended to check the water tank at least quarterly. If contamination is found, the tank should be emptied, cleaned and refilled with new water.

Consumables

To ensure proper functioning of the sterilizer, it is recommended that the wearing parts be replaced according to the schedule below. And periodic inspection of the individual components of the sterilizer in accordance with the following guidelines.

Name	Replacement frequency
HEPA bacteriological filter	Every 1000 cycles or every 12 months
Connection / drain hose	If damage is observed or once a year
Plugs for water / condensate containers	If damage is observed
Door with gasket	Every 1000 cycles or every 12 months

Element to be controlled	Frequency of inspections
Front gasket	weekly or in the event of incorrect operation - performed by the user
Bacteriological filter	every week - performed by the user
Connection / drain hose	weekly or in the event of incorrect operation - performed by the user
Container stoppers	weekly - performed by the user

Data archiving

The course of each performed sterilization is automatically saved on a data carrier (USB stick).

- The memory slot is located on the rear panel of the device.
- It is recommended to periodically archive the data on another medium e.g. computer, laptop.
- The USB flash drive should not be removed from the slot during the process.

Data Viewer software

DataViewer software allows you to view and archive sterilization programs on your computer and print them out.

Installing the software

- The software is delivered with the device and is on a portable drive - USB flash drive.
- To install the software, double click on the software installation file. After this operation an installation window will be displayed regarding language selection.
- After confirming, you must accept the license terms for the software being installed.
- You will then be prompted to place a shortcut to the software on your computer desktop.
- After clicking “Next>”, the folder where the program should be installed is selected.
- After making your selection, click “Next”. By clicking the Install button, the Data Viewer software will be installed.


Program structure and main functionalities

The main window consists of three main areas.

Open	Process no. 00186	Date 07-05-21	Start 12:52:46	Synchronise
Search	Program 134 Fast (S)		Stop 13:00:22	Print/PDF Report

File Search

- ST01-CH-21-10835
 - 2021-04-22
 - 2021-04-23
 - 2021-04-26
 - 2021-04-27
 - 2021-04-28
 - 2021-04-29
 - 2021-04-30
 - 2021-05-04
 - 2021-05-05
 - 2021-05-06
 - 2021-05-07
 - LOG00176.dat
 - LOG00177.dat
 - LOG00178.dat
 - LOG00179.dat
 - LOG00180.dat
 - LOG00181.dat
 - LOG00182.dat
 - LOG00183.dat
 - LOG00184.dat
 - LOG00185.dat
 - LOG00186.dat**
 - LOG00187.dat
 - LOG00188.dat
 - 2021-05-10
 - 2021-05-11
 - 2021-05-12
 - 2021-05-13
 - 2021-05-14
 - 2021-05-17

Stage	Time [mm:ss]	Step [mm:ss]	p [bar]	t [°C]	t calc. [°C]	
0	Locking door					
1	Heating up	00:01	00:01	-0.067	89.9	89.9
2	Vacuum I	00:03	01:15	-0.866	94.6	94.6
3	Steam Pulse I					
4	Vacuum II					
5	Steam Pulse II					
6	Vacuum III					
7	Pressurizing	01:19	01:05	2.203	135.9	135.9
8	Sterilizing	02:25	03:40	2.201	136.1	136.5
9	Depressurizing	06:06	01:24	-0.398	113.9	113.9
10	Drying					
11	Equalizing	07:31	00:03	-0.122	114.4	114.4
12	Unlocking door	07:35		-0.077	114.6	114.6

Sterilization summary	STERILIZATION TEMPERATURE,	134	Time	03:40
Max Temperature	136,2	Max Pressure	2,24	
Min Temperature	135,8	Min Pressure	2,19	

Sterilization Complete

Note:

Save

Tree of all processes that have been synchronized with the flash drive - have been sorted by execution dates

Temperature and pressure graph with main autoclave and process data (date and number).

Data on the duration and achieved parameters of each process step.

The most important sterilization parameters.

Possibility to save notes for each process

The function buttons, eg “PDF report”, which will enable printing the protocol from the process, are marked in dark blue.

Drop-down menu:

By clicking on the File window, we have access to the following options:

- Loading the saved process flow from a flash drive or from another location.
- Print a saved process.
- Create a report to a PDF file.
- Exporting data to a file in order to send a package with data to the manufacturer in case of problems.
- Data export to CSV format.

By clicking on the Tools window, we have access to the following options:

- Synchronization of all files with saved processes after selecting memory location on pendrive (this process, depending on the number of sterilizations and tests performed can take up to several seconds).
- Searching for any saved process from the database.
- Adding your own logo to PDF reports.

Search

The program allows you to search for processes by:

- Date range.
- Process number.
- Process type.
- Process outcome (success/fail).

PDF report

The program allows you to generate a report on each process performed by the autoclave. It contains all the necessary process data and the steriliza-

tion result.

Error codes

The table below contains error codes that may appear when working with the sterilizer

No	Error code	Description	Recommendations
1	“Chamber over temperature”	Maximum chamber temperature exceeded	Contact with the service
2	“Steam gen. Over temperature”	Steam generator temperature too high	<ul style="list-style-type: none"> • Too high weight of sterilized instruments - repeat the process with less instruments (max. 0.5 kg) • Contact with the service
3	“Process over temperature”	Process temperature too high	Contact with the service
4	“Overpressure error”	Pressure error	Contact with the service
5	“Sterilization pressure too low”	Drying pressure too low	<ul style="list-style-type: none"> • Check that there is water in the bottle with the blue hose • Correct the position of the blue hose so that the end is completely submerged in water. Add a sinker to eliminate the problem in the future • Check that the water supply hose (blue) is not damaged (After correcting the position / replacing the hose or refilling with water, restart the machine) • Contact with the service
6	“Sterilization temp. Too low”	Sterilization temperature too low	<ul style="list-style-type: none"> • Check the water level in the bottle with the blue hose • Check that the red tube is not pointing up along its entire length, creating the so-called air trap • Contact with the service
7	“Too high pressure during drying”	Drying pressure too high	<ul style="list-style-type: none"> • Make sure the red drain hose is not immersed in water. The hose must not be kinked, the liquid must flow down by gravity • Check that the weight of the sterilized instruments is not too high • Contact with the service
8	“Too many steam pulses / no water”	Too many steam pulses. No feed water.	<ul style="list-style-type: none"> • Check the water connection to the “water in” connector • Check the distilled water level in the feed water tank (blue plug) • Check that the weight of the load does not exceed the allowable weight • Contact with the service
9	“Drainage error”	Clogged drain	<ul style="list-style-type: none"> • Clogged drain Check the level of the wastewater and the connection of the hoses • Check the level of the used water in the bottle with the red cap. If the bottle is full, discard the used water • Check that the red hose is not kinked and that it is pointing downwards along its entire length • Check that there is no debris in the outlet opening (inside the chamber) • Contact with the service
10	“Chamber heating error”	Chamber heating error	<ul style="list-style-type: none"> • Mains voltage too low - consult an electrician for the location where the autoclave is to be installed • Contact with the service

11	“Steam generator heating error”	Steam generator error	<ul style="list-style-type: none"> • Repeat the process • Contact with the service
12	“Prevacuum fail / check outlet”	Vacuum pump / drain fault	<ul style="list-style-type: none"> • Check the level of used water in the bottle with the red cap and pour out the excess • Check that the red hose is not submerged or kinked • The red hose must point downwards along its entire length, no section may point upwards • Clean the chamber door seal • Contact with the service
13	“Power failure”	Momentary voltage loss during operation	<ul style="list-style-type: none"> • Restart the device and make sure it is properly plugged into the outlet • Consult an electrician competent for the site where the autoclave is installed to check the installation
14	“Pressure during standby”	Overpressure while waiting	<ul style="list-style-type: none"> • Restart your device • Contact with the service
15	“Locking door error”	Door lock error	Contact with the service
16	“Unlocking door error”	Door unlock error	<ul style="list-style-type: none"> • Turn off the autoclave and turn it on again, start the process and stop it after a few seconds. There must be no overpressure in the chamber, i.e. the following information must be displayed: “READY / Chamber is safe to open” • Contact with the service
17	“Valve V3 / HEPA filter error”	V3 valve / HEPA filter error	<ul style="list-style-type: none"> • Replace the HEPA filter • Contact with the service
18	“Pressure sensor error”	Pressure sensor error	Contact with the service
19	“USB disc error / Change Disc”	Writing error on pendrive - damage to medium	Copy the contents from your current flash drive - buy and install a new one
20	“Min. Chamber temperature”	Chamber temperature too low during the process	Contact with the service
21	“Chamber temperature sensor failure”	Chamber temperature sensor failure	Contact with the service
22	“Steam gen. Temp. Sensor failure”	Steam generator temperature sensor failure	Contact with the service
23	“Process temp. Sensor failure”	Process temperature sensor failure	Contact with the service
24	“Autoclave has too low”	Autoclave temperature too low / temperature sensor error	<ul style="list-style-type: none"> • Leave the device switched off for 3 hours at room temperature • Contact with the service
31	“Internal flash error”	Internal memory error	Contact with the service

Messages

“Aborted by user”	Process interrupted by the user. Non-sterile cartridge if interrupted during or before the sterilization process	This message appears when the user terminates the process. This does not mean that there is a malfunction. Start a new process.
“Vacuum test failed”	Vacuum test error	Contact with the service
“No USB memory”	No USB memory stick	Check the USB port and mount the memory. Contact with the service.
“Equalizing pressure”	Pressure when stationary Equalization of pressure to atmospheric	<ul style="list-style-type: none"> • The message occurs in certain cases as a result of natural processes • If the message appears frequently, contact the service center

“Overpressure during standby”	Hypertension in standby mode	The reason for this error is that the hot sterilizer is left with the chamber closed (e.g. overnight). As the sterilizer cools down a vacuum is created in the chamber which causes a startup error
		Wait until the device has equalized the pressure automatically - the message will disappear automatically

TECHNICAL DATA

Technical data	Maxima Turbo
Power supply	110-120 V/60Hz
Installed power	1.6 kW max.
Maximum electric current consumption	15 A
Operating pressure	2.1 bar / 30.5 psi
Maximum pressure	2.45 bar max / 35.53 psi
Maximum process temperature	137°C (278°F)
Process chamber capacity	2.7 l / 0.7 gal
Mass	15 kg / 33 lb
Process chamber dimensions (LxWxH)	292 x 192 x 45 mm / 11.4 x 7.5 x 1.8 in
External device dimensions (LxWxH)	561 x 252 x 162 mm / 22 x 9.8 x 6.3 in
Protection rating	IP20
Noise level	49dB(A)
Process data archiving	USB drive

SERVICE

For items that require service, including warranty service, please send to Henry Schein® Handpiece & Small Equipment Solutions.

Henry Schein Handpiece & Small Equipment Solutions
620 S. Placentia Ave.
Placentia, CA 92870 USA
Phone # 800-235-6337 Prompt 1

WARRANTY

This product is designed for use in a dental office and this warranty is not applicable to other uses. This product is warranted against manufacturing defects in materials and workmanship for 2 years from date of purchase.

Henry Schein® will repair or replace product at its sole discretion. This warranty shall be limited to replacement or repair of the product or its parts and shall not extend to any other claims, including but not limited to loss of profit, cost of removal or replacement, incidental or consequential damages or other similar claims arising from use of this product. Damage to product resulting from acts of God, faulty installation, misuse, tampering, accident, abuse, negligence, alterations, unauthorized repairs or problems unrelated to materials and workmanship are not covered by this warranty.

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