

AUTOCLAVES & STERILIZERS

MAXIMA TURBO CLASS B STERILIZER

INDUSTRY-LEADING SPEED

Triple vacuum cycle sterilization and dry time in 15 minutes.





*Compared to a standard autoclave

For more information, see product Instructions for Use.

Order today at HENRYSCHEINDENTAL.COM/HSBRAND









AUTOCLAVES & STERILIZERS

MAXIMA TURBO CLASS B STERILIZER

Class B Sterilizer is defined by a pre-sterilization vacuum cycle. The class B is considered to be the highest class of autoclave and can be used to sterilize all loads.

Advantages

Industry-leading Speed

Triple-vacuum technology leaves your instruments perfectly sterile and dry in 15 minutes

Designed for Efficiency

Fits up to eight pouches per 15-minute cycle

Quiet Operation

Performs equivalent to a quiet conversation, providing no disruption at your practice

Save Space in Style

Units are portable, stackable, and can conveniently fit into small spaces

Emergency Sterilization Mode

Sterilize unwrapped instruments for immediate use in only seven minutes

Features

- · Easy-to-read, high-resolution digital display
- · User-friendly one-click, one-button navigation
- · Simple installation process
- · Built-in USB port for transferring your cycle data right to your PC
- · Uses distilled water from external tank



MARKET-LEADING 2-YEAR WARRANTY!

INCLUDED IN BOX: Maxima Turbo Sterilizer | Water and Condensate Connection Cables | Quick Start Guide | USB Drive | Instructions for Use HEPA Filter

UNWRAPPED FAST / 273°F

WRAPPED VACUUM / 273°F

WRAPPED VACUUM / 250°F

Capacity	
Max. packages amount (circa)	8 pouches
Pouch dimensions	3.5" x 9"
Dimensions	
Unit size (inches)	9.8"L x 6.3"H x 22"D

Also Available

Maxima PowerClean 210 Ultrasonic Cleaner

(570-2241)



Maxima Turbo Class B Sterilizer

(570-3362)

Related Product

Water Jug 1 Gallon, (2) per pack

(570 - 3376)

Order today at HENRYSCHEINDENTAL.COM/HSBRAND





